


teknova:

DESCRIPTION: 200L Bioprocess Bag



NOTE: 200L CHAMBER DESIGNED TO FIT A
57 GALLON POLY DRUM WITH LEVER LOCK LID

teknova:

DESCRIPTION: 200L Bioprocess Bag

BILL OF MATERIALS

NO.	PART#	QTY	UDM	DESCRIPTION
30	970310	1	EA	CONNECTOR, STRAIGHT, 3/8 X 3/8, PP
29	120040	1	EA	PORT FLANGE, 3/8 OPPOSING HB, ENGAGE 8402
28	160025	2	EA	PKG MTRL, 4" X 6" POLYBAG, TAMPER-EVIDENT
27	940131	6	EA	CABLE TIE, SUB MINIATURE, HTSN
26	940130	28	EA	CABLE TIE, 5-1/2" STANDARD, HTSN
25	960068	3	EA	PINCH CLAMP, STD SMALL, 1/4 OD, PP
24	960070	3	EA	PINCH CLAMP, STD MEDIUM, 1/2 OD, PP
23	960067	3	EA	PINCH CLAMP, STD LARGE, 5/8 OD, PP
22	920359	1	EA	COVER, VINYL, Ø0.562 X L1"
21	110116	3	EA	FITTING, FEMALE LUER CAP, W/O STEM, PC
20	110342	3	EA	FITTING, NEEDLE-FREE INJECTION, FLL X 1/8 HB, PC
19	970758	1	EA	CONNECTOR, CROSS, 1/4 ALL, PP
18	970585	1	EA	CONNECTOR, REDUCING TEE, 3/8 X 3/8 X 1/4, PP
17	970114	1	EA	CONNECTOR, REDUCER, 1/2 > 3/8, PP
16	110006	1	EA	FITTING, COLDER MPC CAP FOR INSERT, PC
15	110005	1	EA	FITTING, COLDER MPC INSERT, 3/8 HB, PC
14	110007	1	EA	FITTING, COLDER MPC PLUG FOR BODY, PC
13	110003	1	EA	FITTING, COLDER MPC BODY, 3/8 HB, PC
12	100484	1	EA	FILTER, SARTORIUS, 0.2/0.1µm PES, 5445358M8G-00
11	190064	0.3	FT	TUBING, ADVANTAFLEX TPE, 1/2 X 3/4
10	190011	0.8	FT	TUBING, ADVANTAFLEX TPE, 1/8 X 1/4
9	190017	2.8	FT	TUBING, ADVANTAFLEX TPE, 1/4 X 7/16
8	190019	15.5	FT	TUBING, ADVANTAFLEX TPE, 3/8 X 5/8
7	150153	3	EA	PORT TUBE ASSY, Ø3/8 X 1/8 HB, FX# & PP
6	150154	3	EA	PORT TUBE ASSY, Ø3/8 X 1/4 HB, FX# & PP
5	950140	2	EA	CHAMBER, FLEXTREME, 150mL, FlexClr™ BPM
4	950142	1	EA	CHAMBER, FLEXTREME, 500mL, FlexClr™ BPM
3	120044	2	EA	PORT FLANGE, 3/8 HOSE BARB, ENGAGE 8402
2	990235	24.0	FT	FILM, FLEXCLR BP (OUTER), W22-1/2
1	990288	24.0	FT	FILM, FLEXCLR M (INNER), W22-1/2

teknova:

DESCRIPTION: 200L Bioprocess Bag