




NOTES	
A	Refer to subsequent print for BOM



teknova:

DESCRIPTION: 50L x 4 Manifold



BILL OF MATERIALS			
ITEM	DESCRIPTION	QTY.	MATERIAL
1	50 L, 4-Port, 3D Biocontainer	4	Proprietary
2	Protective Cover, Medium	1	Vinyl
3	Meissner Capsule Filter CL2STW0.1-222	1	PP/Proprietary
4	3/8" x 5/8" Tubing, 12" Length	2	Proprietary
5	3/8" x 3/8" x 1/4" HB T-Fitting	1	Polypropylene
6	1/4" x 7/16" Tubing, 4" Length	4	Proprietary
7	1/4" HB X-Fitting	1	Polypropylene
8	1/4" x 7/16" Tubing, 10" Length	3	Proprietary
9	Ratchet Pinch Clamp, Medium	10	Polypropylene
10	Male Luer Cap	7	Polycarbonate
11	1/4" HB Swabable Valve / Needleless Injection Site	7	Polycarbonate/Silicone
12	250 mL End-Ported Biocontainer	2	Proprietary
13	500 mL End-Ported Biocontainer	1	Proprietary
14	3/8" HB T-Fitting	3	Polypropylene
15	Ratchet Pinch Clamp, Large	12	PBT
16	3/8" x 5/8" Tubing, 48" Length	14	Proprietary
17	Coupling Body, 3/8" HB, MPC	4	Polycarbonate
18	Sealing Plug, MPC	4	Polycarbonate
19	1/4" x 7/16" Tubing, 6" Length	4	Proprietary
20	Coupling Insert, 3/8" HB, MPC	4	Polycarbonate
21	Sealing Cap, MPC	4	Polycarbonate
22	Tamper Evident Re-closable Poly-bag (not pictured)	15	Polyethylene
23	Bubble Wrap Poly-Bag, Medium (not pictured)	1	Polyethylene
24	Bubble Wrap Poly-Bag, Large (not pictured)	4	Polyethylene
25	Cable Tie, Standard (not pictured)	55	Nylon
26	Cable Tie, Small (not pictured)	15	Nylon
27	Outer Packaging, Poly-Bag (not pictured)	2	Polyethylene

NOTES	
B	All connections secured with a single cable tie (except at 'C')
C	Connection identified secured with dual cable ties, located 180° opposite each other
D	Filter BOM item no. 3 to be placed in a bubble wrap poly-bag, which is then sealed
E	Secure poly-bags over swabable valve / needleless injection site & cap and secure with single cable ties
F	Secure poly-bags over MPC fittings and secure with single cable ties
G	Tubing BOM items no. 16 to be coiled, placed together (per set attached to each biocontainer) in bubble wrap poly-bags, which are then sealed

teknova:

DESCRIPTION: 50L x 4 Manifold