

Custom Pre-filled Bioprocess Bags

Single-use bags are ideal for streamlining R&D or clinical large-volume and scale-up applications. Choose stock solutions, buffers, and media in various bag formats, including pillow bags, rocking bioreactor bags, and 3D process tank bags. Alternatively, have your custom formulations delivered in bioprocess bags.

- Select from 2D and 3D bags (5–200 L) or totes from 500 L
- Minimize the risk of contamination
- Improve reproducibility
- Reduce bottle waste



buffers and solutions • animal-free media • dry culture media • liquid culture media • pre-poured agar plates

Common inventoried items available in bioprocess bags

Description	Bag Size	Part No.
LB Broth	20 L	L8020
Terrific Broth Plus MOPS	20 L	T8130
1X PBS	200 L	P2020
1X PBS, Endotoxin Tested	20 L	2P1825
Ethanol 20%, Endotoxin Tested	20 L	2E7525
Sodium Hydroxide (0.1 N NaOH) and Sodium Chloride (1 M NaCl)	20 L	2S8915
Water, DNase/RNase/Endotoxin Tested	20 L 200 L	2W1230 2W0200

For a complete list of inventoried items or to enquire about custom bioprocess bags, contact us at info@teknova.com.

teknova.com

For Research Use Only.

© 2021 Alpha Teknova, Inc. All rights reserved. Teknova is a trademark of Alpha Teknova, Inc., and is registered in the USA. All other marks are the property of their respective owners. FL10014 08/2021

TEKNOVA
science matters